



Recommended Devices and Operating System for PointClickCare Skin and Wound Product

Operating system	V11 or higher Please upgrade your device before installing PointClickCare's Skin and Wound application.
Recommended Devices	iPhone - 6S 6S Plus 7 7 Plus 8 8 Plus X XR XS XS Max 11 11 Pro iPad - Pro 2 nd Gen Pro 3 rd Gen
Supported Devices Devices with a flash are preferred.	iPhone - 6 6 Plus SE iPod Touch - 6 th Gen 7 th Gen iPad - Regular 6 th Gen Air 2 Mini 4 Mini 5

Web Dashboard Specification

Internet Explorer	V11 or higher**		
Edge	Most stable recent version**		
Google Chrome	Most stable recent version**		
Mozilla Firefox	Most stable recent version**		
Apple Safari	Most stable recent version**		
TLS	V1.2		
Whitelisted Domains	On Firewall for HTTPS 443: *.swiftmedical.io https://swiftmedical-prod.s3.amazonaws.com/* https://swiftmedical-prod-background.s3.amazonaws.com/* https://pcc-swiftmedical-io-background.s3.amazonaws.com/* https://s3.amazonaws.com/static.dashboard.swiftmedical.io/* wss://ws.pusherapp.com https://api-mt1.pusher.com/ https://js.pusher.com/ https://api.mixpanel.com		

IP addresses for PointClickCare.com domain

The following table lists the US & Canada PointClickCare Public IP Addresses:

USA	Internap Chicago, USA	Internap Chicago, USA	Cogent Chicago, USA
Public IP Block	64.74.10216/28	64.74.109.128/25	38.103.114.112/29
IP Range	64.74.102.16 - 64.74.102.31	64.74.109.128 - 64.74.109.255	38.103.114.112 - 38.103.114.119
	Reliance Chicago	Reliance Chicago	
Public IP Block	67.200.121.64/26	67.200.121.16/29	
IP Range	67.200.121.64 - 67.200.121.127	67.200.121.16 - 67.200.121.23	

CANADA	Bell Toronto, Canada	Bell Toronto, Canada	Bell Toronto, Canada
Public IP Block	72.28.70.96/27	72.28.70.128/28	72.28.70.224/28
IP Range	72.28.70.96 - 72.28.70.127	72.28.70.143 - 72.28.70.143	72.28.70.224 - 72.28.70.239